


File : L:\QC Data\GCMS-06 Minimal\DATA\121319\12-.D  
Operator :  
Acquired : 13 Dec 19 09:51 pm using AcqMethod CHICK\_2019.M  
Instrument : GCMS 6  
Sample Name : ALR-069S-CN  
Misc Info :  
Vial Number : 147

 **AccuStandard®**  
Leader in Analytical Reference Standards  
Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 80c, 17c/min to 340c, 8min  
GC Parameters: Cons. Split, 12psi constant flow  
Split 20:1, 1uL inj.; GC/MS; INJ 270c

