

File :Q:\GCMS-08 Minimal\DATA\20191031\15-D  
Operator : KM / JP / ZA  
Acquired : 31 Oct 2019 18:37 using AcqMethod KERLA-2017.M  
Instrument : GCMS 8  
Sample Name: ALR-113S  
Misc Info : Ethyl Paraben  
Vial Number: 134



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

