

File :Q:\GCMS-08 Minimal\DATA\20230828\34- [JP].D
Operator : Organic QC Lab
Acquired : 28 Aug 23 09:45 pm using AcqMethod Hy_PAH_24cm_S200.M
Instrument : GCMS 8
Sample Name: AOCs-005N
Misc Info : in CS2
Vial Number: 23



Column: DB-5, 30m, 0.25ID, 0.25um
Oven: 80c, 2m, 30c/m to 320, final 2min
GC para: Pulse split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj, GCMS, Inj 280c

