

File :Q:\GCMS-08 Minimal\DATA\20210917\17- YX.D  
Operator : Organic QC Lab  
Acquired : 17 Sep 2021 16:22 using AcqMethod Hy\_PAH\_24cm\_S100.M  
Instrument : GCMS 8  
Sample Name:APP-9-004  
Misc Info : Aectophenone in DCM  
Vial Number: 103



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

