

File :D:\MassHunter\GCMS\1\DATA\012723\07 - (sak).D
Operator : Organic QC Lab
Acquired : 27 Jan 2023 09:50 using AcqMethod CHICK_2019_13MIN.M
Instrument : GCMS 6
Sample Name : APP-9-066
Misc Info :
Vial Number: 113



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

