

File : Q:\GCMS-04\DATA\112818\80-  
Operator : RFXM  
Acquired : 29 Nov 2018 15:13 using AcqMethod VOC-100SPLIT DECKER.M  
Instrument : GCMS4  
Sample Name: APP--9-211  
Misc Info : RT  
Vial Number: 60

Column: DB-624 UI, 30m x 0.25mm ID x 1.4 $\mu$ m  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1 $\mu$ l inj.; GC/MS; Injector temp 240°C

