

File :Q:\GCMS-08 Minimal\DATA\20220920\18- -JP-.D
Operator : Organic QC Lab
Acquired : 20 Sep 22 02:34 pm using AcqMethod Hy_PAH_24cm_S100.M
Instrument : GCMS 8
Sample Name: AS-E0079
Misc Info : Phenanthrene; in meoh
Vial Number: 4



Column: DB-5, 30m, 0.25 ID, 0.25um
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

