

File : Q:\GCMS-07\DATA\10241801\35-  
Operator : K.D  
Acquired : 25 Oct 2018 00:27 using AcqMethod KERLA-2Split-2uL.M  
Instrument : 5977B MSD  
Sample Name: AS-E0107  
Misc Info :  
Vial Number: 69

Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

