

File : B:\QC Data\GCMS-04 Minimal\DATA\100719\19-.D
Operator : ZA / RM / KM
Acquired : 7 Oct 2019 20:14 using AcqMethod VOC-200SPLIT DECKER.M
Instrument : GCMS4
Sample Name : AS-E0175
Misc Info :
Vial Number: 61



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

