

File :Q:\GCMS-04 Minimal\DATA\010423A\89- RM.D
Operator :
Acquired : 7 Jan 2023 20:02 using AcqMethod VOC-100 SPLIT #H2 28MIN(HOT).M
Instrument : GCMS4
Sample Name: AS-E0337
Misc Info :
Vial Number: 22



Column: DB-624UI, 20m x 0.18mm ID x 1.0 μ m
Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
GC Parameters: 1 μ l inj.; GC-MS; Injector temp 250 $^{\circ}$ C carrier gas Hydrogen

