

File : Q:\GCMS-08 Minimal\DATA\20230714\45- [JP].D
Operator : Organic QC Lab
Acquired : 15 Jul 23 10:22 am using AcqMethod Hy_PAH_24cm_S150.M
Instrument : GCMS 8
Sample Name: AS-E0411
Misc Info : Acetophenone; in meoh
Vial Number: 39



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

