

File :D:\MassHunter\GCMS\1\data\20231003\23- (JP).D
Operator : Organic QC Lab
Acquired : 03 Oct 23 04:35 pm using AcqMethod Hy_PAH_24cm_S200.M
Instrument : GCMS 8
Sample Name: AS-E0458
Misc Info : N-Nitrosopiperidine; in meoh
Vial Number: 14



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

