

File : Q:\GCMS-04 Minimal\DATA\120821\28- RM.D
Operator :
Acquired : 9 Dec 2021 1:56 using AcqMethod VOC-100SPLIT .M
Instrument : GCMS4
Sample Name: AS-E1104
Misc Info : RT
Vial Number: 14

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

