

File :Q:\GCMS-04 Minimal\DATA\040220\39-D
Operator : RFXM
Acquired : 3 Apr 2020 2:03 using AcqMethod VOC-200SPLIT HOT280.M
Instrument : GCMS4
Sample Name: ASTM-P-0050
Misc Info :
Vial Number: 69



Column: ZB-624Plus, 30m, 0.25 ID, 1.40um
Oven Parameters: 35c, 5m, 11c/min to 60c, 22c/min to 280c, 2min
GC Parameters: Split 100:1, 1uL inj.; GC/MS INJ 240c

