

File :K:\QC Data\GCMS-04 Minimal\DATA\052821\15-RM.D
Operator :
Acquired : 28 May 2021 19:01 using AcqMethod VOC-100SPLIT DECKER.M
Instrument : GCMS4
Sample Name: ASTM-P-0091-10-10X
Misc Info :
Vial Number: 33

Column: DB-624 UI, 30m x 0.25mm ID x 1.4 μ m
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1 μ l inj.; GC/MS; Injector temp 240°C

