

File : B:\QC Data\GCMS-04 Minimal\DATA\091219\77-.D
Operator : RFXM/KM
Acquired : 13 Sep 19 10:04 pm using AcqMethod VOC-SCAN 24 DECKER.M
Instrument : GCMS4
Sample Name: BF-MEOH-5X
Misc Info
Vial Number: 33



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

