

File :D:\MassHunter\GCMS\1\data\20220425\06 - SAK.D
Operator : Organic QC Lab
Acquired : 25 Apr 2022 13:26 using AcqMethod Hy_PAH_24cm_S100.M
Instrument : GCMS 8
Sample Name : BPA-G-S
Misc Info :
Vial Number: 71



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

