

File :D:\MassHunter\GCMS\1\data\20201230\15- JP.D
Operator : Organic QC Lab
Acquired : 30 Dec 20 03:19 pm using AcqMethod Hy_Pest_15cm_S10.M
Instrument : GCMS 8
Sample Name: C-023S-TP
Misc Info : 2,3,5-trichlorobiphenyl in isooctane
Vial Number: 62



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

