

File :Q:\GCMS-06 Minimal\DATA\012521\13- -DC.D
Operator : Organic QC Lab
Acquired : 25 Jan 2021 17:27 using AcqMethod CHICK_2019_13MIN.M
Instrument : GCMS 6
Sample Name: C-030S
Misc Info : 2,4,6-Trichlorobiphenyl in Isooctane
Vial Number: 37



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

