

File :M:\GCMS-06 Minimal\DATA\102121\07- ALR.D
Operator : Organic QC Lab
Acquired : 21 Oct 2021 18:42 using AcqMethod CHICK_2019_13MIN.M
Instrument : GCMS 6
Sample Name: C-035S-TP
Misc Info : 3,3',4-Trichlorobiphenyl in Isooctane
Vial Number: 3



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

