

File :D:\MassHunter\GCMS\1\data\20210826\14- (JP).D
Operator : Organic QC Lab
Acquired : 27 Aug 21 01:22 am using AcqMethod Hy_Pest_15cm_S5.M
Instrument : GCMS 8
Sample Name: C-036S-TP
Misc Info : 33'5-Tri-PCB in iso
Vial Number: 64



Column: DB-5, 30m, 0.25ID, 0.25um
Oven: 80c, 2m, 30c/m to 320, final 2min
GC para: Pulse split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj, GCMS, Inj 280c

