

File :D:\MassHunter\GCMS\1\data\053019\37-.D
Operator : D.C
Acquired : 31 May 19 06:26 am using AcqMethod CHICK_2019.M
Instrument : GCMS 6
Sample Name: C-044S-TP
Misc Info : C-044 in Isooctane
Vial Number: 48



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

