


File : L:\QC Data\GCMS-06 Minimal\DATA\110719\45-.D
Operator : K M/YX/ZA
Acquired : 08 Nov 2019 06:31 using AcqMethod CHICK_2019_S100.M
Instrument : GCMS 6
Sample Name : C-068N
Misc Info : RT
Vial Number : 146

 **AccuStandard**
Leader in Analytical Reference Standards
Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

