

File :Q:\GCMS-08 Minimal\DATA\20210129\40- (JP).D
Operator : Organic QC Lab
Acquired : 30 Jan 2021 00:50 using AcqMethod Hy_Pest_24cm_S2_2ulinj.M
Instrument : GCMS 8
Sample Name: C-169S-TP
Misc Info : 3,3',4,4',5,5'-HexaPCB; in isooctane
Vial Number: 86



Column: DB-5, 30m, 0.25ID, 0.25um
Oven: 80c, 2m, 30c/m to 320, final 2min
GC para: Pulse split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj, GCMS, Inj 280c

