

File : Q:\GCMS-07 Minimal\DATA\10242201\12- JP-.D
Operator : Organic QC department
Acquired : 24 Oct 2022 11:22 using AcqMethod KERLA-2017_Short_15Min.M
Instrument : 5977B MSD
Sample Name: C-203S-TP
Misc Info : 22'344'5'6-Octa-PCB; in iso
Vial Number: 141



Column: DB-5, 30m, 0.25 ID, 0.25um
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

