

Chrom Perfect Chromatogram Report



Sample Name: C-216-260-WL-5X-5ML

File: C:\CPData\DATA\20220520\222051334 JP.0000.RAW

Date (start): 5/20/2022 5:09:35 PM

Column: 110c. 0min, 15c/min to 325, Final 2 min

Parameters: Splitless; 2ul; ECD: 325c; Inj: 250c

notes:


Method: Ni_FAST_ECD.MET Setpoint File: Ni_FAST_ECD.M68

Calibration File:

Inst: GC #7

Oper: JJP

Vial = 2



| | |
|-------|------------------------------------|
| 8.38 | Tetrachloro-meta-xylene (877-09-8) |
| 10.36 | Aroclor 1016 (12674-11-2) |
| 14.32 | Aroclor 1260 (11096-82-5) |
| 16.84 | Decachlorobiphenyl (2051-24-3) |