



AccuStandard®

Sample Name: C-221-ST-2 ; 500 ug/g; dil 30x in iso

File: C:\CPData\DATA\20211019\221101319 jp.0000.RAW

Date (start): 10/19/2021 11:46:52 AM

Column: 110c, 0min, 15c/min to 325, Final 2 min

Parameters: Splitless; 2ul; ECD: 325c; Inj: 250c

notes:

Method: Ni_FAST_ECD.MET

Setpoint File: Ni_FAST_ECD.M68

Calibration File:

Inst: GC #7 Oper: JJP

Vial = 2

— C:\CPData\DATA\20211019\221101319 jp.0000.RAW; 500 ng/ml in iso

