



AccuStandard®

Sample Name: C-242-ST-2(Aroclor-1242 500ug/g trsfmr oll; dil 50x in iso

File: C:\CPData\DATA\20220829\222081380 JP(1).0000.RAW

Date (start): 8/29/2022 4:00:17 PM

Column: 80c,1m,10c/m to 320,2m

Parameters: Inj 320c, ECD 320c

notes:

Method: Ni_PCB_S50_ECD.MET

Setpoint File: Ni_PCB_S50_ECD.M68

Calibration File:

Inst: GC #7

Oper: JJP

Vial = 4

Response - Millivolts

