

File :Q:\GCMS-06 Minimal\DATA\120419\06-.D  
Operator :  
Acquired : 04 Dec 2019 17:42 using AcqMethod CHICK\_2019\_2S\_2uL.M  
Instrument : GCMS 6  
Sample Name: C-248S-TP  
Misc Info : Aroclor 1248 in Isooctane  
Vial Number: 14



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

