

File :D:\MassHunter\GCMS\1\DATA\050422\11 -SAK.D
Operator : Organic QC Lab
Acquired : 04 May 2022 09:20 using AcqMethod CHICK_2019_2S_2uL.M
Instrument : GCMS 6
Sample Name: C-260S-M
Misc Info :
Vial Number: 49



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

