

File :Q:\GCMS-06 Minimal\DATA\052019\32-.D  
 Operator : D.C  
 Acquired : 21 May 19 02:30 am using AcqMethod CHICK\_2017\_S2\_2uL.M  
 Instrument : GCMS 6  
 Sample Name: C-CS-04  
 Misc Info : PCB Congener Mix #4 in Isooctane  
 Vial Number: 106



Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
 Oven Program: 80c, 17c/min to 340c, 8min  
 GC Parameters: Cons. Split, 12psi constant flow  
 Split 2:1, 2uL inj.; GC/MS; INJ 270c

