

File :Q:\GCMS-08 Minimal\DATA\20220715\30- =JP=.D
 Operator : Organic QC Lab
 Acquired : 15 Jul 2022 21:21 using AcqMethod Hy_PCB_2.5cm_S2.M
 Instrument : GCMS 8
 Sample Name: C-WNN
 Misc Info : PCBs in iso
 Vial Number: 135



Column: DB-5, 30m, 0.25ID, 0.25um
 Oven: 80c, 2m, 30c/m to 320, final 2min
 GC para: Pulse split 40psi 0.5m, 10psi constant flow
 Split 20:1, 1ul inj, GCMS, Inj 280c

