

File : Q:\GCMS-04 Minimal\DATA\093021\28- RM.D

Operator :

Acquired : 1 Oct 2021 2:35 using AcqMethod VOC-100SPLIT DECKER.M

Instrument : GCMS4

Sample Name: CLP-004-10X

Misc Info :

Vial Number: 46

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

