

File :Q:\GCMS-04\DATA\101118\19-D
 Operator : RFXM
 Acquired : 11 Oct 2018 19:55 using AcqMethod VOC-100SPLIT DECKER.M
 Instrument : GCMS4
 Sample Name: CLP-022-R2
 Misc Info :
 Vial Number: 87



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

