

File :Q:\GCMS-08 Minimal\DATA\20240115\18- (JP).D
Operator : Organic QC Lab
Acquired : 15 Jan 24 06:06 pm using AcqMethod Hy_PAH_24cm_S200.M
Instrument : GCMS 8
Sample Name: CLP-HC-A-R
Misc Info : in dcm
Vial Number: 11



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

