

File : M:\GCMS-07 Minimal\DATA\06222201\13-  
 Operator : Organic QC department  
 Acquired : 22 Jun 2022 11:15 using AcqMethod KERLA-2017.M  
 Instrument : 5977B MSD  
 Sample Name: CP-ORE-03-R1  
 Misc Info :  
 Vial Number: 1

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
 Split 20:1, 1µl inj.; GC/MS; Injector 270°C

