

File :D:\MassHunter\GCMS\1\DATA\08012201\02-KD.D  
 Operator : Organic QC department  
 Acquired : 01 Aug 22 07:54 am using AcqMethod KERLA-VOC-SV.M  
 Instrument : 5977B MSD  
 Sample Name: CP-TER-MIX-003  
 Misc Info :  
 Vial Number: 1

Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
 Oven Program: 80c, 17c/min to 340c, 8min  
 GC Parameters: Cons. Split, 12psi constant flow  
 Split 20:1, 1uL inj.; GC/MS; INJ 270c

