

File :B:\QC Data\GCMS-04 Minimal\DATA\081919\21-.D
Operator : RFXM
Acquired : 19 Aug 2019 22:19 using AcqMethod VOC-200SPLIT HOT.M
Instrument : GCMS4
Sample Name: D-3606R1-IS2-07
Misc Info :
Vial Number: 67



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

