

File :Q:\GCMS-06\DATA\032519\06a-.D  
Operator : D.C.  
Acquired : 25 Mar 19 05:23 pm using AcqMethod M502-Split200.M  
Instrument : GCMS 6  
Sample Name: D-4815-03  
Misc Info : ASTM D-4815 Oxygenates  
Vial Number: 73



Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 35c, 25c/min to 250c, 2min  
GC Parameters: Cons. Split, 120mL/min constant flow  
Split 200:1, 1uL inj.; GC/MS; INJ 200c

