

File :L:\QC Data\GCMS-04 Minimal\DATA\030521B\07--RM.D
Operator : ORGANIC QC DEPT.
Acquired : 5 Mar 2021 16:09 using AcqMethod VOC-100SPLIT DECKER.M
Instrument : GCMS4
Sample Name: D-4815-VT
Misc Info :
Vial Number: 12



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

