

File : Q:\GCMS-04 Minimal\DATA\110222\42- RM.D
Operator :
Acquired : 3 Nov 2022 12:42 using AcqMethod VOC24-100SPLIT#H2 25MIN 255DEG.M
Instrument : GCMS4
Sample Name: D-5501-12-01
Misc Info :
Vial Number: 30



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm
Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

