

File : Q:\GCMS-04 Minimal\DATA\111319A\16-.D
Operator : RM / KM / ZA
Acquired : 13 Nov 19 4:50 pm using AcqMethod VOC-SCAN 24 DECKER.M
Instrument : GCMS4
Sample Name: D-5501-12-02
Misc Info : dil. 10:1 in IPA
Vial Number: 56



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

