

File :K:\QC Data\GCMS-04 Minimal\DATA\033020\15-.D
Operator : RFXM
Acquired : 30 Mar 2020 21:54 using AcqMethod VOC-SCAN 24 200SPLIT.M
Instrument : GCMS4
Sample Name: D-5599-QC-IS-2ML
Misc Info :
Vial Number: 33



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

