

File : Q:\GCMS-04\DATA\120418\61-
 Operator : RFXM/KM
 Acquired : 5 Dec 2018 17:17 using AcqMethod VOC-100SPLIT DECKER.M
 Instrument : GCMS4
 Sample Name: D-5769-ADD-IS-R
 Misc Info :
 Vial Number: 81



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

