

File :Q:\GCMS-06\DATA\032919\26-.D
 Operator : D.C.
 Acquired : 30 Mar 2019 06:00 using AcqMethod CHICK_2017_VOC-SV_S200.M
 Instrument : GCMS 6
 Sample Name: D-5769-CAL-05
 Misc Info :
 Vial Number: 16



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
 Split 20:1, 1µl inj.; GC/MS; Injector 270°C

