

File : C:\msdchem\1\DATA100118A\05-
 Operator : RFXM
 Acquired : 1 Oct 2018 16:55 using AcqMethod VOC-100SPLI
 Instrument : GCMS4
 Sample Name: D-5769-QC-15ML
 Misc Info :
 Vial Number: 58

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

