

File : Q:\GCMS-04 Minimal\DATA\081922\22- RM.D  
 Operator :  
 Acquired : 20 Aug 2022 2:52 using AcqMethod VOC-100SPLIT#H2 20MIN175DEG.M  
 Instrument : GCMS4  
 Sample Name: DIN38407-9-BENZ  
 Misc Info : RT  
 Vial Number: 53

Column: DB-624UI, 20m x 0.18mm ID x 1.0µm  
 Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold  
 GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

