

File :Q:\GCMS-08 Minimal\DATA\20220830\20- (JP).D
Operator : Organic QC Lab
Acquired : 30 Aug 2022 14:43 using AcqMethod Hy_PAH_40c_24cm_S2_2ulinj.M
Instrument : GCMS 8
Sample Name: DINEN-12673
Misc Info :in ethanol
Vial Number: 2



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

