

File : Q:\GCMS-08 Minimal\DATA\20191010\BSB\K08-D
Operator : KM / JP
Acquired : 10 Oct 2019 19:58 using AcqMethod C5-C55_S10.M
Instrument : GCMS 8
Sample Name: DRH-002S
Misc Info : D-2887 Calibration Standard in Carbon disulfide
Vial Number: 43



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

