

File :D:\MassHunter\GCMS\1\DATA\122822\11 - 2sak.D
Operator : Organic QC Lab
Acquired : 28 Dec 2022 11:00 using AcqMethod CHICK_2019_VOC-SV_S100.M
Instrument : GCMS 6
Sample Name: DRH-TX-002-D-40X
Misc Info :
Vial Number: 124



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

